

All^{ele}ustrious pmWasabi-Actin

(Cytoskeleton)

Catalog Number: ABP-FP-WACT100

Size: 10ug Price: \$349.00

Introduction

All^{ele}ustrious pmWasabi-actin is a mammalian expression vector that encodes a fusion of mWasabi N-terminal to human beta actin ACTB. It can be a great tool for visualizing cytoskeleton (microfilaments).

All^{ele}ustrious mWasabi is a monomeric green fluorescent protein that can be easily detected using standard GFP filter sets. mWasabi may be used as a direct replacement for EGFP or other GFPs for superior performance, and may be co-imaged with blue and red fluorescent labels without substantial bleed-through.

One of the most conserved proteins, actin can be found in microfilaments, thin filaments and other contractile apparatus in all eukaryotic species. Actin has several isoforms, which are globular proteins of about 43 kD. Actins can be present in extremely high concentrations, e.g. over 100 uM. Actin filaments are linked to alpha-actinin and to membrane through vinculin. The head domain of vinculin associates with E-cadherin via alpha-, beta-, and gamma-catenins. The tail domain of vinculin binds to membrane lipids and to actin filaments (see Alleleustrious fusions with these proteins).

Source

Engineered variant of mTFP1, originally derived from *Clavularia sp.* coral.

Recommended Use

mWasabi has been optimized for use with standard GFP/FITC filter sets.

Features

- About 2-fold brighter than EGFP
- Similar photostability to EGFP
- Uses standard filter sets
- Can be co-imaged with blue and red FPs or dyes
- Mammalian expression vector ready to transfect your favorite cells
- Low sensitivity to acidic pH (fluorescence pKa=4.3)
- True monomer that will not aggregate or cause nonspecific interactions

Reconstitution

10 µg provided in lyophilized powder form. Reconstitute with 10 µL of nuclease-free water for a final concentration of 1 µg/µL.

Storage

Store at -20°C or at -80°C for long-term preservation.

Human CMV Immediate-Early

Promoter (CMV Promoter).....1-595

mWasabi.....613-1320

Actin.....1321-2469

SV40 PolyA Signal.....2629-2679

bla Promoter.....3218-3322

Ampicillin Resistance Gene.....3307-4167


pUC Origin.....4316-4958

Upstream Sequencing Primer:

Universal CMV Promoter Primer

Downstream Sequencing Primer:

SV40 Primer:GCTTT ATTG TGAA TTTGT GATGC TATTG C



References: Ai H, Olenych SG, Wong P, Davidson MW, Campbell RE. Hue-shifted monomeric variants *Clavularia* cyan fluorescent protein: identification of the molecular determinants of color and applications in fluorescence imaging. *BMC Biology*. 2008 Mar; 6:13. Shaner NC, Patterson GH, Davidson MW. Advances in fluorescent protein technology. *J Cell Sci*. 2007 Dec 15;120(Pt 24):4247-60. Ai HW, Hazelwood KL, Davidson MW, Campbell RE. Fluorescent protein FRET pairs for ratiometric imaging of dual biosensors. *Nature Methods*. 2008 5(5): 401-03. Ai HW, Henderson JN, Remington SJ, Campbell RE. Directed evolution of a monomeric, bright, and photostable version of *Clavularia* cyan fluorescent protein: structural characterization and applications in fluorescence imaging. *Biochem J*. 2006. Shaner NC, Steinbach PA, Tsien RY. A guide to choosing fluorescent proteins. *Nat Methods*. 2005 2(12):905-09. Holmes KC, Popp D, Gebhard W, Kabsch W. Atomic model of the actin filament. *Nature* 1990, 347, 21-2. Berg T; S Mousavi, Malerod L, Kjeker R. "Clathrin-dependent endocytosis". *Biochem J* 1990 377: 1-16. Khaitlina SY Functional specificity of actin isoforms. *Int Rev Cytol* 2001 202:35-98.

SEQUENCE

```

1 TAGTATTAA TAGTAATCAA TTACGGGGTC ATTAGTTCAT AGCCCATATA TGGAGTTCG
61 CGTTACATAA CTTACGGTAA ATGGCCCCC TGGCTGACCG CCCAACGACC CCCGCCATT
121 GACGTCAATA ATGACGTATG TTCCCATAGT AACGCCAATA GGGACTTTC ATTGACGTCA
181 ATGGGTGGAG TATTTACGGT AAAGTCCCA CTGGCAGTA CATCAAGTGT ATCATATGCC
241 AAGTACGCC CCTATTGACG TCAATGACGG TAAATGGCCC GCCTGGCATT ATGCCAGTA
301 CATGACCTTA TGGGACTTTC CTACTTGGCA GTACATCTAC GTATTAGTCA TCGCTATTAC
361 CATGGTGATG CGGTTTGGC AGTACATCAA TGGGCGTGA TAGCGGTTG ACTCACGGGG
421 ATTTCCAAGT CTCCACCCCA TTGACGTCAA TGGGAGTTG TTTTGGCACC AAAATCAACG
481 GGACTTTCCA AATGTTCGTA ACAACTCCG CCAATTGACG CAAATGGGG GTAGCCGTGT
541 ACGGTGGGAG GTCTATATAA GCAGAGCTGG TTTAGTGAAC CGTCAGATCC GCTAGCGCTA
601 CCGTTCGCCA CCATGGTGAG CAAGGGCGAG GAGACCACA TGGGCGTAAT CAAGCCGAC
661 ATGAAGATCA AGCTGAAGAT GGAGGGCAAC GTGAATGGCC ACGCCTTCGT GATCGAGGG
721 GAGGGCGAGG GCAAGCCCTA CGACGGCAC AACACCATCA ACCTGGAGGT GAAGGAGGGA
781 GCCCCCTGC CTTTCTCCTA CGACATTCTG ACCACCGCGT TCAGTTACGG CAACAGGGCC
841 TTCACCAAGT ACCCCGACGA CATCCCAAC TACTTCAAGC AGTCTTCCC CGAGGGCTAC
901 TCTTGGGAGC GCACCATGAC CTTTCGAGGAC AAGGGCATCG TGAAGGTGAA GTCCGACATC
961 TCCATGGAGG AGGACTCCTT CATCTACGAG ATACACCTCA AGGGCGAGAA CTTCCCCCC
1021 AACGGCCCCG TGATGCAGAA GGAGACCACC GGCTGGGACG CCTCCACCGA GAGGATGTAC
1081 GTGCGCGACG GCGTGTGTA GGCGGACGTC AAGATGAAGC TGCTGTGGA GGGCGCGGG
1141 CACCACCGCG TTGACTTCAA GACCATCTAC AGGGCCAAGA AGGCGGTGAA GCTGCCCGAC
1201 TATCATTTTG TGGAACCGC CATCGAGATC CTGAACCACG ACAAGGACTA CAACAAGGTG
1261 ACCGTTTACG AGATCGCCGT GGCCCGCAAC TCCACCGACG GCATGGACGA GCTGTACAAG
1321 TCCGGACTCA GATCTCGAGC TATGGATGAT GATATCGCCG CGCTCGTCGT CGACAACGGC
1381 TCCGGCATGT GCAAGGCCGG CTTTCGCGGGC GACGATGCCC CCCGGGCCGT CTTCCCTCC
1441 ACCGAGCGGG GCCCAGGCA CACAGGCGTG ATGGTGGGCA TGGGTGAGAA GGATCCTAT
1501 GTGGGCGACG AGGCCAGAG CAAGAGAGGC ATCCTCACC TGAAGTACC CATCGAGCAC
1561 GGCATCGTCA CCAACTGGGA CGACATGGAG AAAATCTGGC ACCACACCTT CTACAATGAG
1621 CTGCGTGTGG CTCCCGAGGA GCACCCCGTG CTGCTGACCG AGGCCCCCTT GAACCCCAAG
1681 GCCAACCGGG AGAAGATGAC CCGATCATG TTTGAGACCT TCAACACCCC AGCCATGTAC
1741 GTTGCTATCC AGGCTGTGCT ATCCCTGTAC GCCTCTGGCC GTACCACTGG CATCGTGATG
1801 GACTCCGGTG ACGGGGTAC CCACACTGTG CCCATCTACG AGGGGTATGC CCTCCCCAT
1861 GCCATCTGCG GTCTGGACCT GGCTGGCCGG GACCTGACTG ACTACCTCAT GAAGATCCTC
1921 ACCGAGCGGG GCTACAGCTT CACCACCACG GCCGAGCGGG AAATCGTGCG TGACATTAAG
1981 GAGAAGCTGT GCTACGTGCG CCTGGACTTC GAGCAAGAGA TGGCCACGGC TGCTTCCAGC
2041 TCCTCCCTGG AGAAGAGCTA CGAGCTGCCT GACGGCCAGG TCATCACCAT TGGCAATGAG
2101 CGGTTCCGCT GCCCTGAGGC ACTCTTCCAG CTTTCTTCC TGGGCATGGA GTCCTGTGGC
2161 ATCCACGAAA CTACCTTCAA CTCCATCATG AAGTGTGACG TGGACATCCG CAAAGACCTG
2221 TACGCCAACA CAGTGTGTC TGGCGGCACC ACCATGTACC CTGGCATTGC CGACAGGATG
2281 CAGAAGGAGA TCACTGCCCT GGCAACCAGC ACAATGAAGA TCAAGATCAT TGCTCTCCT
2341 GAGCGCAAGT ACTCCGTGTG GATCGGGCGG TCCATCTCTG CCTCGCTGTC CACCTTCCAG
2401 CGATGTGGAG TCAGCAAGCA GGAGTATGAC GAGTCCGGCC CCTCCATCGT CCACCGGAAA
2461 TGCTTCTAGG GATCCACCGG ATCTAGATAA CTGATCATAA TCAGCCATAC CACATTTGTA
2521 AAGTTTTTAC TTGCTTTAAA AAACCTCCA CACCTCCCC TGAACCTGAA ACATAAAATG
2581 AATGCAATTG TTGTTGTTAA CTTGTTTATT GCAGCTTATA ATGGTTACAA ATAAAGCAAT
2641 AGCATCACAA ATTFCACAAA TAAAGCATTT TTTTCACTGC ATCTAGTTG TGGTTTGTCC
2701 AAACATCACTA ATGTATCTTA ACGCGTCTG TAGCGGCGCA TTAAGCGCGG CGGGTGTGGT
2761 GGTACGCGC AGCGTGACCG CTACACTTGC CAGCGCCCTA GCGCCCGCTC CTTTCGCTTT
2821 CTTCCCTTCC TTTCTCGCCA CGTTTCGCGG CTTTCCCGT CAAGCTCTAA ATCGGGGGCT
2881 CCTTGTAGG TCCGATTTA GTGCTTTACG GCACCTCGAC CCAAAAAAAC TFGATTAGGG
2941 TGATGGTTCA CGTAGTGGGC CATCGCCCTG ATAGACGGTT TTTGCGCCCTT TGACGTTGGA
3001 GTCCACGTTT TTTAATAGTG GACTCTTGT CCAAACCTGA ACAACACTCA ACCCTATCTC
3061 GGTCTATTCT TTTGATTTAT AAGGGATTTT GCCGATTTG GCCTATTGGT TAAAAAATGA
3121 CGTATTTTAA CAAAAATTA ACGCGAATTT TAACAAAATA TTAACGCTTA CAATTTAGTT
3181 GGCACCTTTC GGGGAAATGT GCGCGGAACC CCTATTTGTT TATTTTCTA AATACATTCA
3241 AATATGTATC CGCTCATGAG ACAATAACCC TGATAAATGC TCAATAATA TTGAAAAAGG
3301 AAGAGTATGA GTATTCAACA TTTCCGTGTC GCCCTTATTC CTTTTTTTGC GGCATTTTGC
3361 CTTCTGCTTT TGCTCACCC AGAAAACGCTG GTGAAAAGTAA AAGATGCTGA AGATCAGTTG
3421 GGTGCACGAG TGGGTTACAT CGAACTGGAT CTCAACAGCG GTAAGATCCT TGAGAGTTTT
3481 CGCCCGAAG AACGTTTTCC AATGATGAGC ACTTTTAAAG TCTGTCTATG TGGCGCGGTA
3541 TTATCCCGTA TTGACGCGCG GCAAGAGCAA CTCGGTTCGC GCATACACTA TTCTCAGAAT
3601 GACTTGGTTG AGTACTCAC AGTCACAGAA AAGCATCTTA CGGATGGCAT GACAGTAAGA
3661 GAATTATGCA GTGCTGCCAT AACCATGAGT GATAACACTG CGGCCAATT ACTTCTGACA
3721 ACGATCGGAG GACCGAAGGA GCTAACCGCT TTTTTCGACA ACATGGGGGA TCATGTAACT
3781 CGCCTTGATC GTTGGGAACC GGAGCTGAAT GAAGCCATAC CAAACGACGA GCGTGACACC
3841 ACGATGCTG TAGCAATGGC AACCAAGTTG CGCAAACTAT TAACCTGGCA TAACCTACT
3901 CTAGCTTCCC GGCAACAATT AATAGACTGG ATGGAGGCGG ATAAAGTTGC AGGACCACTT
3961 CTGCGCTCCG CCCTTCCGGC TGGCTGGTTT ATGCTGATA AATCTGGAGC CGGTGAGCGT
4021 GGTCTCGCG GTATCATTGC AGCACTGGGG CCAGATGGTA AGCCCTCCC TATCGTAGTT
4081 ATCTACAGGA CGGGGAGTCA GGCAACTATG GATGAACGAA ATAGACAGAT CGTCAAGATA
4141 GGTGCCTCAC TGATTAAGCA TTGGTAACTG TCAGACCAAG TTTACTCATA TATACTTTAG
4201 ATTGATTTAA AACTTCATTT TTAATTTAAA AGGATCTAGG TGAAGATCCT TTTTGATAAT
4261 CTCATGACCA AAATCCCTTA ACGTGAGTTT TCGTTCCACT GAGCGTCAGA CCCCGTAGAA
4321 AAGATCAAAG GATCTTCTTG AGATCCTTTT TTTCTGCGC TAATCTGCTG CTTGCCAACA
4381 AAAAAACCAC CGCTACCAGC GGTGGTTTGT TTGCCGATC AAGAGCTACC AACTCTTTTT
4441 CCGAAGGTAA CTGGCTTCA GAGAGCGCAG ATACCAATA CTGTTCTTCT AGTGTAGCCG
4501 TAGTTAGGCC ACCACTTCAA GAACTCTGTA GCACCGCCTA CATACCTCGC TCTGCTAATC

```



4561 CTGTTACCAG TGGCTGCTGC CAGTGGCGAT AAGTCGTGTC TTACCGGGTT GGA CTCAAGA
4621 CGATAGTTAC CGGATAAGGC GCAGCGGTCG GGCTGAACGG GGGGTTTCGTG CACACAGCCC
4681 AGCTTGGAGC GAACGACCTA CACCGAACTG AGATACCTAC AGCGTGAGCT ATGAGAAAGC
4741 GCCACGCTTC CCGAAGGGAG AAAGGCGGAC AGGTATCCGG TAAGCGGCAG GGTCGGAACA
4801 GGAGAGCGCA CGAGGGAGCT TCCAGGGGGA AACGCCTGGT ATCTTTATAG TCCTGTCGGG
4861 TTTCGCCACC TCTGACTTGA GCGTCGATTT TTGTGATGCT CGTCAGGGGG GCGGAGCCTA
4921 TGGAAAAACG CCAGCAACGC GGCCTTTTTA CGGTTCCCTGG CCTTTTGCTG GCCTTTTGCT
4981 CACATGTTCT TTCCTGCGTT ATCCCCTGAT TCTGTGGATA ACCGTATTAC CGCCATGCAT



Enzyme	Cuts	Positions	Sequence					
AatII	5	125 178 261 447 1109	gacgt/c	BstSFI	6	1167 1932 2728 3848 4526 4717	c/tryag	
Acc113I	2	2350 3613	agt/act	BstSNI	1	340	tac/gta	
Acc16I	1	3871	tgc/gca	BstX2I	12	585 1226 1330 1913 2470 2479	r/gatcy	
AccB1I	9	465 745 1448 1457 1598 2245 2301	g/gyrcc	BstXI	1	762	ccannnn/ntgg	
		2910 4141		BstYI	12	585 1226 1330 1913 2470 2479	r/gatcy	
AccB7I	1	1052	ccannnn/ntgg			3447 3464 4232 4244 4330 4341		
AccBSI	4	1959 2103 2810 3255	gagcgg	BstZI	1	1949	c/ggccg	
AccI	1	1370	gt/mkac	Bsu36I	1	2114	cc/tnagg	
AccIII	1	1321	t/ccgga	Cfr10I	5	600 1048 1396 2846 4009	r/ccggy	
AccsI	3	2649 3134 3145	r/aatty	Cfr9I	1	1421	c/ccggg	
AcyI	8	122 175 258 444 1058 1106 1449	gr/cgyc	CfrI	8	696 1137 1776 1884 1949 2021	y/ggccc	
		3554				2073 3701		
AfeI	1	596	agc/gct	CvnI	1	2114	cc/tnagg	
AflIII	2	2721 4982	a/crygt	DraI	4	2537 3516 4208 4227	ttt/aaa	
AgeI	1	600	a/ccggt	DraII	3	836 1662 1890	rg/gnccy	
AhdI	1	4094	gacnnn/ngtc	DraIII	4	1209 1647 1827 2954	caacnn/gtg	
Alw21I	5	1559 1643 3426 3511 4672	gwgw/c	DrdI	3	2077 2998 4880	gacnnn/ngtcc	
Alw44I	2	3422 4668	g/tgca	DsaI	6	360 611 962 1277 1946 2024	c/crygg	
AlwNI	1	4573	cagnnn/ctg	DsaI	8	696 1137 1776 1884 1949 2021	y/ggccc	
Ama87I	4	889 1334 1421 1633	c/ycgrg			2073 3701		
AocI	1	2114	cc/tnagg	EagI	1	1949	c/ggccg	
Aor51HI	1	596	agc/gct	Eam1104I	3	2057 2127 3305	ctcttc	
ApaLI	2	3422 4668	g/tgca	Eam1105I	1	4094	gacnnn/ngtcc	
ApoI	3	2649 3134 3145	r/aatty	EaRI	3	2057 2127 3305	ctcttc	
AseI	2	7 3919	at/taat	EclHKI	1	4094	gacnnn/ngtcc	
AsnI	2	7 3919	at/taat	EclXI	1	1949	c/ggccg	
Asp700I	1	3494	gaann/nnttc	Eco105I	1	340	tac/gta	
AspEI	1	4094	gacnnn/ngtcc	Eco130I	4	360 611 962 1676	c/cwggg	
AspHI	5	1559 1643 3426 3511 4672	gwgw/c	Eco24I	2	783 2880	grgcy/c	
AspI	1	1813	gacn/ngtcc	Eco255I	2	2350 3613	agt/act	
AtsI	1	1813	gacn/ngtcc	Eco31I	4	636 1047 1719 4027	ggcttc	
AvaI	4	889 1334 1421 1633	c/ycgrg	Eco32I	1	1353	gat/atc	
AviII	1	3871	tgc/gca	Eco47III	1	596	agc/gct	
BalI	2	698 2023	tgg/cca	Eco52I	1	1949	c/ggccg	
BamHI	1	2470	g/gatcc	Eco57I	5	678 1101 1545 3412 4460	ctgaag	
BanI	9	465 745 1448 1457 1598 2245 2301	g/gyrcc	Eco64I	9	465 745 1448 1457 1598 2245 2301	g/gyrcc	
		2910 4141				2910 4141		
BanII	2	783 2880	grgcy/c	Eco81I	1	2114	cc/tnagg	
BbeI	1	1452	ggcgc/c	Eco88I	4	889 1334 1421 1633	c/ycgrg	
BbiII	8	122 175 258 444 1058 1106 1449	gr/cgyc	Eco91I	2	1257 1815	g/gtnacc	
		3554		EcoO109I	3	836 1662 1890	rg/gnccy	
BbsI	1	1434	gaagac	EcoO65I	2	1257 1815	g/gtnacc	
Bbv12I	5	1559 1643 3426 3511 4672	gwgw/c	EcoRV	1	1353	gat/atc	
Bbv16II	1	1434	gaagac	EcoT14I	4	360 611 962 1676	c/cwggg	
BcgI	2	1850 3578	cgannnnntg	EcoT22I	1	5039	atgca/t	
BclI	1	2492	t/gatca	EheI	1	1450	ggc/gcc	
BcoI	4	889 1334 1421 1633	c/ycgrg	ErhI	4	360 611 962 1676	c/cwggg	
BglI	6	90 212 283 1180 1423 3976	gcnnnn/nggc	Esp1396I	1	1052	ccannnn/ntgg	
BglII	1	1330	a/gatct	FauNDI	1	234	ca/tatg	
BpiI	1	1434	gaagac	FbaI	1	2492	t/gatca	
BpmI	4	768 1131 2052 4009	ctggag	FriOI	2	783 2880	grgcy/c	
BpuAI	1	1434	gaagac	FspI	1	3871	tgc/gca	
BsaAI	3	340 1080 2951	yac/gtr	GsuI	4	768 1131 2052 4009	ctggag	
BsaBI	3	2323 2407 2497	gatnn/nnttc	HaeII	5	598 1452 2796 2804 4742	rgcgc/y	
BsaHI	8	122 175 258 444 1058 1106 1449	gr/cgyc	HinII	8	122 175 258 444 1058 1106 1449	gr/cgyc	
		3554				3554		
BsaI	4	636 1047 1719 4027	ggcttc	HincII	3	1152 1371 2598	gty/rac	
BsaMI	2	2585 2684	gaatgc	HindII	3	1152 1371 2598	gty/rac	
BsaOI	5	605 1952 3576 3725 4648	cgry/cg	HpaI	1	2598	gtt/aac	
BsaWI	7	600 1321 1804 2476 3798 4629	w/ccggw	Hsp92I	8	122 175 258 444 1058 1106 1449	gr/cgyc	
		4776				3554		
Bse118I	5	600 1048 1396 2846 4009	r/ccggy	KasI	1	1448	g/gcgc	
Bse21I	1	2114	cc/tnagg	Kpn2I	1	1321	t/ccgga	
Bse8I	3	2323 2407 2497	gatnn/nnttc	Ksp22I	1	2492	t/gatca	
BseAI	1	1321	t/ccgga	Ksp632I	3	2057 2127 3305	ctcttc	
BseRI	5	633 972 1641 2045 2338	gaggg	MamI	3	2323 2407 2497	gatnn/nnttc	
Bsh1285I	5	605 1952 3576 3725 4648	cgry/cg	MfeI	1	2585	c/aattg	
Bsh1365I	3	2323 2407 2497	gatnn/nnttc	MflI	12	585 1226 1330 1913 2470 2479	r/gatcy	
BshNI	9	465 745 1448 1457 1598 2245 2301	g/gyrcc			3447 3464 4232 4244 4330 4341		
		2910 4141		MluI	1	2721	a/cgctg	
BsiEI	5	605 1952 3576 3725 4648	cgry/cg	MluNI	2	698 2023	tgg/cca	
BsiHKAI	5	1559 1643 3426 3511 4672	gwgw/c	Mph1103I	1	5039	atgca/t	
BsiI	2	3430 4814	ctcgtg	MroI	1	1321	t/ccgga	
BsiMI	1	1321	t/ccgga	MroNI	2	1396 2846	g/ccggc	
BsmI	2	2585 2684	gaatgc	MscI	2	698 2023	tgg/cca	
BsoBI	4	889 1334 1421 1633	c/ycgrg	MslI	8	365 613 961 1795 2251 3323 3682	caynn/nrtg	
Bspl3I	1	1321	t/ccgga			3841		
Bspl407I	1	1313	t/gtaca	Msp17I	8	122 175 258 444 1058 1106 1449	gr/cgyc	
Bspl43II	5	598 1452 2796 2804 4742	rgcgc/y			3554		
Bspl9I	3	360 611 962	c/catgg	MspA1I	4	2108 3458 4399 4644	cmg/ckg	
BspCI	1	3725	cgat/cg	MunI	1	2585	c/aattg	
BspEI	1	1321	t/ccgga	Mva1269I	2	2585 2684	gaatgc	
BspHI	4	1907 2186 3254 4262	t/catga	NaeI	2	1398 2848	gcc/ggc	
BsplU11I	1	4982	a/catgt	NarI	1	1449	gg/cgcc	
BsrBI	4	1959 2103 2810 3255	gagcgg	NcoI	3	360 611 962	c/catgg	
BsrBRI	3	2323 2407 2497	gatnn/nnttc	NdeI	1	234	ca/tatg	
BsrDI	5	2098 2270 2333 3858 4040	gcaatg	NgOAI	2	1396 2846	g/ccggc	
BsrFI	5	600 1048 1396 2846 4009	r/ccggy	NgOMI	2	1396 2846	g/ccggc	
BsrGI	1	1313	t/gtaca	NheI	1	591	g/ctagc	
BssAI	5	600 1048 1396 2846 4009	r/ccggy	NsiI	1	5039	atgca/t	
BssSI	2	3430 4814	ctcgtg	NspBII	4	2108 3458 4399 4644	cmg/ckg	
BssTII	4	360 611 962 1676	c/cwggg	NspI	2	1389 4986	rcatg/y	
BstD102I	4	1959 2103 2810 3255	gagcgg	PaeR7I	1	1334	c/tcgag	
BstDSI	6	360 611 962 1277 1946 2024	c/crygg	Pf1MI	1	1052	ccannnn/ntgg	
BstEII	2	1257 1815	g/gtnacc	PinAI	1	600	a/ccggt	
BstH2I	5	598 1452 2796 2804 4742	rgcgc/y	Ple19I	1	3725	cgat/cg	
BstI	1	2470	g/gatcc	Ppu10I	1	5035	a/tgcat	
BstMCI	5	605 1952 3576 3725 4648	cgry/cg	PpuMI	1	1890	rg/gwccy	
BstPI	2	1257 1815	g/gtnacc	PshBI	2	7 3919	at/taat	



Psp1406I	2	3492 3865	aa/cggt	SspI	2	3159 3289	aat/att
Psp5II	1	1890	rg/gwccy	StyI	4	360 611 962 1676	c/cwggg
PspAI	1	1421	c/ccggg	Tth111I	1	1813	gacn/nngtc
PspALI	1	1423	ccc/ggg	Van91I	1	1052	ccannnn/ntgg
PspEI	2	1257 1815	g/gtnacc	VneI	2	3422 4668	g/tgcac
PstNHI	1	591	g/ctagc	VspI	2	7 3919	at/taat
PvuI	1	3725	cgat/cg	XbaI	1	2482	t/ctaga
RcaI	4	1907 2186 3254 4262	t/catga	XcmI	3	1299 1869 2043	ccannnn/nnntgg
SalI	1	1369	g/tcgac	XhoI	1	1334	c/tcgag
SapI	1	2058	gctcttc	XhoII	12	585 1226 1330 1913 2470 2479	r/gatcy
ScaI	2	2350 3613	agt/act			3447 3464 4232 4244 4330 4341	
SfcI	6	1167 1932 2728 3848 4526 4717	c/tryag	XmaI	1	1421	c/ccggg
Sfr274I	1	1334	c/tcgag	XmaIII	1	1949	c/ggccc
SmaI	1	1423	ccc/ggg	XmnI	1	3494	gaann/nnttc
SnaBI	1	340	tac/gta	Zsp2I	1	5039	atgca/t
SspBI	1	1313	t/gtaca				

The following enzymes do not cut:

AatI, Acc65I, AclNI, AflII, ApaI, AscI, Asp718I, AvrII, BanIII, BbrPI, BbuI, BfrI, BlnI, BlnI, Bpu1102I, Bpu14I, Bsa29I, BscI, BseCI, BsePI, BsgI, BsiWI, BsmBI, Bsp106I, Bsp119I, Bsp120I, Bsp1720I, Bsp68I, BspDI, BspMI, BspTI, BspXI, BssHII, Bst1107I, Bst98I, BstBI, Bsu15I, CciNI, CelII, Cfr42I, ClaI, CpoI, Csp45I, CspI, Ecl136II, Eco147I, Eco72I, EcoCRI, EcoNI, EcoRI, Esp3I, FseI, HindIII, KpnI, KspI, LspI, MspCI, NotI, Nrul, NspV, PacI, PaeI, Pfl23II, PmaCI, Pme55I, PmeI, PmlI, PshAI, Psp124BI, PspLI, PspOMI, PstI, PvuII, RsrII, SacI, SacII, SbfI, SexAI, SfiI, Sfr303I, SfuI, SgfI, SgrAI, SmiI, SpeI, SphI, Spil, SrfI, Sse8387I, SseBI, SstI, SstII, Stul, SunI, SwaI, Vha464I



Related products:

Current Alle^{ele}ustious Fluorescent Protein Family Members:

The founding member is mTFP1.

mTFPG3 is a green FP with 3 amino difference from mTFP1. It has a slightly red-shifted emission spectrum and is 1.5 fold brighter compared to EGFP. While being very bright, mTFPG3 can be photobleached within ~5 sec, about 30 times faster than EGFP, suitable for certain cell-based assays that require a bright FP with very short half-life.

mTFP0.7 is a precursor during the evolution of mTFP1. It has photo-switchable properties like Dronpa that cycles between fluorescent and nonfluorescent states. It may be developed into components in PALM/SIM applications.

Basic Vectors

Three vectors are available: pNCS-mWasabi, pmWasabi-N and pmWasabi-C.

Subcellular Marker Vectors

Twenty six vectors are available.

Vectors in Viral Vectors

All plasmid format vectors in Allele's Phoenix Retroviral vector or HiTiter Lentiviral Vectors.

IMPORTANT NOTICE TO PURCHASER: This product may be subjects of pending U.S. and foreign patents. Allele Biotech grants academic/government research institutions a worldwide, non-exclusive, royalty-free, limited license to use this product for non-commercial life science research use only. Such license specifically prohibits the right to sell or otherwise transfer this product or its derivatives to a third party.

Industry For-Profit Institutions: Institutions that wish to use this product are required to obtain a sub-licensing agreement from Allele Biotech. For more information please visit Fluorescent Protein Sublicense, or contact us via phone or email.