



For Research Use Only. Not for Diagnostic or Therapeutic Use.

Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Allele Biotech is strictly prohibited

Website: www.allelebiotech.com
Call: 1-800-991-RNAi/858-587-6645
(Pacific Time: 9:00AM-5:00PM)
Email: oligo@allelebiotech.com

Box 1 | Basic Info

Cat. No.	ABP-PAB-11680
Animal ID	RB0465-0466
Host	Rabbit
Reactivity	Human
Format	Affinity-purified
Accession number	NM_002609
Amount	100 µg

Alternative Name(s): JTK12, PDGFR, CD140B, PDGFR1, PDGF-R-beta

PDGFRB, Phospho-tyrosine PDGFR [Y579] polyclonal antibody

Platelet-derived growth factor receptor beta (PDGFRB) is a cell surface tyrosine kinase receptor that binds to members of the platelet-derived growth factor family. These growth factors (PDGFs) are mitogens for cells of mesenchymal origin. A translocation between chromosomes 5 and 12, that fuses the PDGFRB gene to the ETS leukemia gene, results in chronic myelomonocytic leukemia.

Buffers

Purified rabbit polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column and eluted out with both high and low pH buffers and neutralized immediately after elution then followed by dialysis against PBS.

Immunogen

KLH conjugated synthetic peptide comprised of amino acids 575 - 586 [DGHE[Y*]IYVD] of the human phospho-tyrosine PDGFR [Y579] (PDGFRB) protein.

Application:

Tested by peptide-specific ELISA (1:1,000).

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C. Avoid repeated freeze-thaw cycles.

References:

1. Kovalenko M, Denner K, Sandstrom J, Persson C, Gross S, Jandt E, Vilella R, Bohmer F, Ostman A: Site-selective dephosphorylation of the platelet-derived growth factor beta-receptor by the receptor-like protein-tyrosine phosphatase DEP-1. *J. Biol. Chem.* 275(21): 16219-16226 (2000). Steer EJ, Cross NC: Myeloproliferative disorders with translocations of chromosome 5q31-35: role of the platelet-derived growth factor receptor Beta. *Acta Haematol.* 107(2): 113-122 (2002). Review.
2. Wu H, Windmiller DA, Wang L, Backer JM: YXXM motifs in the PDGF-beta receptor serve dual roles as phosphoinositide 3-kinase binding motifs and tyrosine-based endocytic sorting signals. *J. Biol. Chem.* 278(42): 40425-40428 (2003).