



USP4, Ubiquitin Specific Protease 4 polyclonal antibody

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Allele Biotech is strictly prohibited

Website: www.allelebiotech.com
Call: 1-800-991-RNAi/858-587-6645
(Pacific Time: 9:00AM~5:00PM)
Email: oligo@allelebiotech.com

Box 1 | Basic Info

Cat. No.	ABP-PAB-31070
Animal ID	N/A
Host	Rabbit
Reactivity	Human
Format	Purified
Accession number	NM_003363
Amount	50µg

Alternative Name(s):

N/A

Peptide affinity purified anti-Human Ubiquitin Specific Protease 4 (USP4) antibody

Buffers

Purified rabbit polyclonal antibody supplied in PBS, pH 7.7 with 0.01% (W/V) sodium azide. This antibody is purified through a peptide immunoaffinity column and neutralized immediately after elution then followed by dialysis against PBS.

Immunogen


KLH conjugated synthetic peptide of internal domain of human USP4.

Application

N/A

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C. Avoid repeated freeze-thaw cycles.



Thymus