



ATP, High quality, standard rATP suitable for large scale IVT.

Description

Ribonucleotide tri-phosphates (rNTP) are ribonucleotides with 3 phosphate groups. rNTPs are the building blocks of RNA synthesis. Allele offers ATP (nucleoside 5'-triphosphate) rNTP solution, which is suitable for use in *in vitro* transcription reactions with SP6, T3, and T7 RNA polymerases for production of mRNA transcripts. ATPs are supplied in convenient, ready-to-use solutions at a concentration of 100 mM.

Quality Assurance:

Each batch of reagents is vigorously tested for consistency and stability, and certified 99% pure and free of endonucleases and modifying agents.

Website: www.allelebiotech.com
Call: 1-800-991-RNAi/858-587-6645
(Pacific Time: 9:00AM~5:00PM)
Email: oligo@allelebiotech.com

For Technical Support:
Email: oligo@allelebiotech.com

Box 1 | Content

Cat. No.	
ABP-PP-NTPSTDA	
Quantity	Concentration
50µl	100mM

Notes: _____

