

# pLanYFP-NT (ABP-FP-YNNCS)

CMV Promoter ..... 1-595

LanYFP ..... 613-1323

Multiple Cloning Sites ..... 1324-1377


SV40 poly(A) Signal ..... 1554-1604

Ampicillin Resistance Gene ..... 2232-3092

pUC Origin ..... 3241-3883

Sequencing Primer:

LanYFP-F(ABP-FP-YFPFPRM) and SV40pA-R



```

661 CCGGTCGCCA CCATGGTGAG CAAGGGCGAG GAGGATAACA TGGCCTCTCT CCCAGCGACA
721 CATGAGTTAC ACATCTTTGG CTCCTTCAAC GGTGTGGACT TTGACATGGT GGGTCGTGGC
      >>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>
1321 AACTTCAAGG AGTGGCAAAA GGCCTTTACC GATGTGATGG GCATGGACGA GCTGTACAAG
1381 AGATCTCGAG TCGGATCCTC GAATTCGGTA CCCTCTAGAC AAGCTTGCGG CCGCTAA

```

## Multiple Cloning Sites